

neodisher® EM

Detergent component – liquid concentrate

Improves the cleaning performance of alkaline and acidic detergents in particular when removing solid particles, paraffins, oils, and greases in aqueous solutions. With emulsifying, dispersing, and foam control properties.

Article	Conts	Order No.
① neodisher® EM	5 l	10 27 63033

neodisher® FA

„Classical“ alkaline cleaner – liquid concentrate – free of surfactants and oxidants

Excellent cleaning properties with high material compatibility. Removes blood, protein, serum, culture medium residues, and gelatin.

Article	Conts	Order No.
② neodisher® FA	5 l	10 27 61133

neodisher® LaboClean A 8

Alkaline universal detergent – powder – free of surfactants

For universal application, excellent cleaning performance and high dirt carrying properties. Removes blood, serum, food residues, pharmaceuticals, dyes, paints, gelatine, and other organic residues.

Article	Conts	Order No.
③ neodisher® LaboClean A 8	10 kg	10 27 61276

neodisher® LaboClean FLA

Alkaline universal detergent – liquid concentrate – free of surfactants, phosphor connections and oxidants.

For universal application. Excellent cleaning performance. Removes stubborn oil-based residues and other organic contaminations.

Article	Conts	Order No.
④ neodisher® LaboClean FLA	10 l	10 27 61230

neodisher® LaboClean FT

Alkaline detergent – liquid concentrate – with oxidation effect, free of surfactants

Specifically for residue-free cleaning of glassware for cell and tissue cultures, pharmacy glass bottles and infusion flasks, utensils for microbiology and nuclear medicine. Removes residues of culture media, blood, proteins, radioactive contamination, as well as other organic residues.

Article	Conts	Order No.
⑤ neodisher® LaboClean FT	12 kg	10 27 61030

neodisher® LaboClean UW

Alkaline detergent – powder – free of phosphates

Especially suitable for the automated cleaning of laboratory glassware and utensils in water analysis laboratories, in serological institutes and in the phosphates industry, where residues of phosphorous compounds would adulterate analysis results.

Article	Conts	Order No.
⑥ neodisher® LaboClean UW	10 kg	10 27 61976

neodisher® LM 3

Alkaline detergent – liquid concentrate – free of phosphates

For manual cleaning and pre-cleaning of laboratory glass and for cleaning pipettes and viscosimeters in special dishwashers. Particularly well suited for cleaning laboratory glasses water, phosphate, and metal laboratories and for use in labs carrying out enzymatic and serological examinations. Also for use in machines. Low-foaming.

Article	Conts	Order No.
⑦ neodisher® LM 3	10 l	10 27 64030

neodisher® MA

Alkaline detergent – powder – free of oxidants

Excellent cleaning performance combined with high material compatibility.

Article	Conts	Order No.
⑧ neodisher® MA	10 kg	10 27 61776

neodisher® MediClean forte

Alkaline universal detergent – liquid concentrate – contains alkalinity donors, enzymes and surfactants

Excellent material compatibility, universally applicable. Residues of dried and denatured blood are reliably removed. A redeposition of protein residues is prevented.

Article	Conts	Order No.
⑨ neodisher® MediClean forte	5 l	10 27 65033



①



②



③



④



⑤



⑥



⑦



⑧



⑨

neodisher® N

Acid neutralising agent and detergent – liquid concentrate – Free of surfactants

For neutralising any carried over alkaline residues on laboratory glassware or for acid pre-cleaning resp. cleaning. Lime or other acid-soluble residues are easily removed. Easy to rinse off. Based on phosphoric acid.

Article	Conts	Order No.
10 neodisher® N	5 l	10 27 62133
11 neodisher® N	25 kg	10 27 62126

neodisher® Z

Acid neutralising agent and detergent – liquid concentrate – Free of surfactants and phosphates

For neutralising any carried over alkaline residues on laboratory glassware and for removing acid-soluble residues in laboratories where residues of phosphates, nitrogen and surfactants would adulterate sensitive analyses and where particular gentleness on material is vital. Based on organic acids.

Article	Conts	Order No.
12 neodisher® Z	5 l	10 27 62233
13 neodisher® Z	10 l	10 27 62230

weigoman®

Hand disinfectant – ready-made solution – bactericidal, tuberculocidal, fungicidal, virucidal to a limited degree, effective against rota and noro viruses

Alcoholic hand disinfectant for topical application to skin. Excellent skin tolerability, with highly effective lipid replenishers. Suitable for sensitive skin. VAH listed, entered in the IHO list of virucidal products. Dosage by means of the triformin® dispenser 1.0 and the triformin® dispenser 0.5.

Article	Conts	Order No.
14 weigoman®	500 ml	10 27 69097
15 weigoman®	1000 ml	10 27 69059

weigoman® perfume-free

Hand disinfectant – ready-made solution – baktericidal, tuberculocidal, fungicidal, virucidal to a limited degree, effective against rota and noro viruses

Unscented alcoholic hand disinfectant for topical application to skin. Excellent skin tolerability, with highly effective lipid replenishers. Suitable for sensitive skin. VAH listed, entered in the IHO list of virucidal products. Dosage by means of the triformin® dispenser 1.0 and the triformin® dispenser 0.5.

Article	Conts	Order No.
16 weigoman®	500 ml	10 27 69197
17 weigoman®	1000 ml	10 27 69159

neoform® Rapid

Bactericidal, yeasticidal, fungicidal, tuberculocidal, limited spectrum virucidal activity, active against enveloped viruses, active against Norovirus. Alcohol-based solution for disinfecting surfaces of medical devices, surfaces in the laboratory area as well as surfaces and equipment in the pharmaceutical and cosmetics area. Application: wiping method. VAH listed and included in the IHO list of virucidal products.

Article	Conts	Order No.
18 neoform® Rapid	750 ml	10 27 64347

neoform® wipes RTF

Disposable non-woven wipes dispenser system for the use of surface disinfectants. Disinfection of surfaces, equipment, medical devices and medical inventory in combination with suitable Dr. Weigert surface disinfectants in hospitals, medical practices, old people's and nursing homes, in chemical and pharmaceutical areas, laboratories as well as in food processing companies and communal kitchens. Also suitable for high-risk areas such as neonatology or haemato-oncology according to VAH recommendation. Versatile, approved pursuant to the 4-field test (EN 16615).

Article	Conts	Order No.
19 neoform® wipes RTF	115 pcs.	10 27 61388

