

Volumetric Flasks, Class B, VITLAB

Volumetric flasks, PMP, class B

Highly transparent. With ring mark individually calibrated to „In“. Class B tolerances according to DIN EN ISO 1042. Thermal stress up to 121°C (autoclaving) does not permanently exceed the tolerance limit. To preserve markings, cleaning at no higher than 60°C is recommended.

Vol. ml	Tol. ± ml	Height mm	Neck	Order No.	PU pc/Pack
10 ¹⁾	0.08	90	GL 18	10 85 31455	6
25 ¹⁾	0.08	115	GL 18	10 85 31465	6
50 ¹⁾	0.12	150	GL 18	10 85 31472	6
① 100 ¹⁾	0.20	180	GL 18	10 85 31477	6
250 ¹⁾	0.30	235	GL 25	10 85 31484	5
500 ¹⁾	0.50	270	GL 25	10 85 31488	4
1000 ¹⁾	0.80	310	GL 32	10 85 31492	3
10 ²⁾	0.08	90	NS 10/19	10 85 31155	6
25 ²⁾	0.08	115	NS 10/19	23 85 11165	6
50 ²⁾	0.12	150	NS 12/21	23 85 11172	6
② 100 ²⁾	0.20	180	NS 14/23	23 85 11177	6
250 ²⁾	0.30	235	NS 19/26	23 85 11184	5
500 ²⁾	0.50	270	NS 19/26	23 85 11188	4
1000 ²⁾	0.80	310	NS 24/29	23 85 11192	3

1) with screw cap made of PP

2) with stopper NS made of PP

Volumetric flasks, PP, class B

Highly transparent. With ring mark individually calibrated to „In“. Class B tolerances according to DIN EN ISO 1042. Thermal stress up to 60°C does not permanently exceed the tolerance limit. To preserve markings, cleaning at no higher than 60°C is recommended.

Vol. ml	Tol. ± ml	Height mm	Neck	Order No.	PU pc/Pack
10 ¹⁾	0.08	90	GL 18	10 85 31255	6
25 ¹⁾	0.08	115	GL 18	10 85 31265	6
50 ¹⁾	0.12	150	GL 18	10 85 31272	6
③ 100 ¹⁾	0.20	180	GL 18	10 85 31277	6
250 ¹⁾	0.30	235	GL 25	10 85 31284	5
500 ¹⁾	0.50	270	GL 25	10 85 31288	4
1000 ¹⁾	0.80	310	GL 32	10 85 31292	3
10 ²⁾	0.08	90	NS 10/19	10 85 31355	6
25 ²⁾	0.08	115	NS 10/19	10 85 31365	6
50 ²⁾	0.12	150	NS 12/21	10 85 31372	6
④ 100 ²⁾	0.20	180	NS 14/23	10 85 31377	6
250 ²⁾	0.30	235	NS 19/26	10 85 31384	5
500 ²⁾	0.50	270	NS 19/26	10 85 31388	4
1000 ²⁾	0.80	310	NS 24/29	10 85 31392	3

1) with screw cap made of PP

2) with stopper NS made of PP

